



ALFAMEDIC s.r.o.
Tr. 5. května č. 458
373 72 LIŠOV
Czech Republic

E-mail: sales@alfamedic.cz
WWW: <http://www.alfamedic.cz>
Tel.: +420 387 994 664
Fax: +420 387 994 603

ELECTRICAL SUCTION PUMP ES-93



The device is suitable for bed-side suction, for use at operating theatres, at operating theatres in maternity hospitals, at ICU – intensive care units for pathologic newborns, at paediatric departments, at observation units in ambulance cars, etc. It can be used in all situations where devices with simple operation, high reliability and bulkiness are necessary. Noiseless and well-proven diaphragm aggregate with incorporated power unit is used in the device.

The whole device is incorporated in plastic housing; 2 stock bottles for evacuated fluid are fastened at the side of housing. Exchangeable bacteriological bed is built in the device. Suction flask can be with clamp fastened at roll-away rack or at guide rail in ambulance car. Power supply of the device is 230V / 50Hz, it is also supplied with accessories for 12 V car power supply.

Set for continuous suction (e.g. pneumothorax suction) is supplied with adjusted head with possibility of continuous regulation of negative pressure using water column.

TECHNICAL PARAMETERS:

Device:	class 1, type B
Max. negative pressure at suction:	-60 kPa
Supply voltage:	230V / 50Hz, 12V
Dimension HxWxL:	250 x 250 x 250 mm
Weight:	approx. 3 kg
Max. input:	45W
Shielding degree:	EN 55 011
Protection:	2x T315mA / 230V