



ALFAMEDIC s.r.o.
Tr. 5. května č. 458
373 72 LIŠOV
Czech Republic

E-mail: sales@alfamedic.cz
WWW: <http://www.alfamedic.cz>
Tel.: +420 387 994 664
Fax: +420 387 994 603

PHOTOTHERAPY FL-400

Phototherapy lamp FL - 400 is used for treating hyperbilirubinemia of the newborn children. The lamp is installed on a mobile stand with the possibility of height adjusting with simple rotation. The illumination body is side adjustable. These adjusting possibilities make possible the treating of patients in classic beds and incubators. The light source is the electric discharge RVIM 400. The discharge is during operation constantly cooled by the ventilator, which extends its life. There is a switch and an operation hours counter placed on the control case.

FUNCTION:

- height and side adjusting of the illuminating body
- switching the electric discharge lamp
- measuring the operation hours of the discharge

CARE:

The device FL - 400 does not require any special care. It is possible to clean with common detergent.

TECHNICAL PARAMETERS:

Measures:	950 x 650 x 1850 mm
Weight:	cca 40 kg
Min./max. light source adjusting:	1100 - 1800 mm
Light source:	discharge RVIM 400 - blue
Voltage:	230V / 50Hz
Max. power consumption:	500 VA
Fuses:	2 x T10 A
Classification:	Class I.

