



ALFAMEDIC
Tr. 5. května č. 458
373 72 LIŠOV
Czech Republic

E-mail: sales@alfamedic.cz
WWW: http://www.alfamedic.cz
Tel.: +420 38 7994664
Fax: +420 38 7994603

PASTEURISER / HEATER MP-05

Suitable for all workplaces where breast milk is prepared and stored. Combination of pasteurisation / heating mode enables to use the device in most of newborn, nursery and paediatric departments. The device operation is simplified to single choice of working mode

- **pasteurisation**
- **heating**

The device is ready for use after inserting basket with bottles into water bath (up to 64pcs standard baby bottles or 36pcs 500ml bottle).

Heating bath temperature regulation and time needed for pasteurisation or heating of breast milk are controlled automatically, in relation to milk temperature of the inserted bottles.

Termination of chosen mode is indicated by acoustic and light signal. Water bath heating is turned off at the same time. Stainless steel container, baby bottles holder and device housing made of steel material guarantee the observance of sanitary code during breast milk processing at all departments.



TECHNICAL PARAMETERS:

Power supply: 230V / 50 Hz
Power input: max 1800W
Weight: cca 55kg, without filling
Heater: 38,5°C
Pasteur.: 62,5°C
Capacity: 64pcs of baby bottles or
36pcs of 500ml bottles
Operation: touch keyboard
Signalling: optically / acoustic
Charge: approx. 35 l of water
Dimension: 700 x 700 x 795 mm

